



COVID-19 & Influenza A+B & RSV & ADV & HMPV Antigen Combo Test Kit (LFA)

CE 3018




Operational Use Video

**For in vitro diagnostic use only. For self-testing.
Please read the instruction carefully before use.**


[Intended use]

This product is used for the qualitative detection of COVID-19, influenza A, influenza B, respiratory syncytial virus (RSV), adenovirus (ADV) and human metapneumovirus (HMPV) antigen in human nasal swab specimens. It is a non-automated rapid test method for infection. This test is authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals. Individuals who test positive should seek follow up care with their physician or healthcare provider as additional testing may be necessary. Users under the age of 15 should complete the test with supervision of an adult.

[Materials and Components]





[Preparation before the test]



[Test Procedure]

Allow test device extraction reagent and specimens to equilibrate to room temperature (15 ~ 30 °C) prior to testing. Please keep the temperature at 15 ~ 30 °C and the humidity at 20%-80% during the whole test.



[Interpretation of test results]

Positive result:

- (1) Flu A Positive: Two red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (A), and the other is in the control zone (C).
- (2) Flu B Positive: Two red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (B), and the other is in the control zone (C).
- (3) COVID-19 Positive: Two red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (COV), and the other is in the control zone (C).
- (4) Flu A/Flu B Positive: Three red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (A), one is in the test zone (B), and the other is in the control zone (C).
- (5) Flu A/COVID-19 Positive: Three red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (A), one is in the test zone (COV), and the other is in the control zone (C).
- (6) Flu B/COVID-19 Positive: Three red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (B), one is in the test zone (COV), and the other is in the control zone (C).

(7) Flu A/Flu B/COVID-19 Positive: Four red bands appear in the Flu A/Flu B/COV detection zone. One is in the test zone (A), one is in the test zone (B), one is in the test zone (COV), and the other is in the control zone (C).

(8) RSV Positive: Two red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (RSV), and the other is in the control zone (C).

(9) ADV Positive: Two red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (ADV), and the other is in the control zone (C).

(10) HMPV Positive: Two red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (HMPV), and the other is in the control zone (C).

(11) RSV/ADV Positive: Three red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (RSV), one is in the test zone (ADV), and the other is in the control zone (C).

(12) RSV/HMPV Positive: Three red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (RSV), one is in the test zone (HMPV), and the other is in the control zone (C).

(13) ADV/HMPV Positive: Three red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (ADV), one is in the test zone (HMPV), and the other is in the control zone (C).

(14) RSV/ADV/HMPV Positive: Four red bands appear in the RSV/ADV/HMPV detection zone. One is in the test zone (RSV), one is in the test zone (ADV), one is in the test zone (HMPV), and the other is in the control zone (C).

[Summary]

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main Manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Once infected with the COVID-19 virus, you may be hospitalized and some complications may occur.

Influenza, usually called "flu", is an acute respiratory infectious disease caused by influenza viruses. It is highly contagious and is spread mainly through coughing and sneezing. It usually breaks out in spring and winter. Divided into influenza A virus, influenza B virus and influenza C virus. Influenza A virus has strong variability, followed by influenza B virus, and influenza C virus is very stable, so influenza A virus is more serious and prevalent than influenza B virus.

Respiratory syncytial virus (RSV) is a common, and very contagious, virus that infects the respiratory tract of most children before their second birthday. Nearly half of all children become infected by RSV in their first year of life. It is also the major viral cause of nosocomial illness in children already hospitalized for other reasons. In children hospitalized with RSV infection, it is believed to be the most common viral cause of death in children younger than 5 years, particularly in children younger than one year. RSV infection can cause cold-like symptoms, including a cough and runny nose, which usually last 1 to 2 weeks. Respiratory syncytial virus spreads through the air, like after a cough or a sneeze, and through direct contact like touching. People are typically infected with RSV for the first time as an infant or toddler and nearly all children are infected before their second birthday. However, repeat infections may occur throughout life, and people of any age can be infected. Infections in healthy children and adults are generally less severe than among infants and older adults with certain medical conditions.

Adenovirus(ADV) is a kind of double-stranded DNA virus without envelope, which mainly causes respiratory tract, eye and gastrointestinal tract infection. Symptoms of respiratory illness caused by Adenovirus infection range from the common cold syndrome to pneumonia, croup and bronchitis. Patients with compromised immune systems are especially susceptible to severe complications of Adenovirus infection. Adenovirus is transmitted by direct contact, fecal-oral transmission and occasionally waterborne transmission. Some types are capable of establishing persistent asymptomatic infections in tonsils, adenoids and intestines of infected hosts and shedding can occur for months or years.

Human metapneumovirus (HMPV) was first reported in 2001 in the Netherlands and since then has been found to be distributed world wide and is associated with a significant proportion of the respiratory infections in early infancy and childhood, but is present in all age groups. It is RNA virus that belongs to the Metapneumovirus genus within the Pneumovirinae subfamily of the Paramyxoviridae family, which can cause respiratory infection. HMPV infection is a typical respiratory infection in infants. HMPV, like other common respiratory viruses including RSV and Influenza, causes seasonal epidemics centred around the winter months, and symptoms are very similar to those experienced by patients infected with RSV, including most commonly fever with upper respiratory symptoms such as cough and nasal congestion. However, HMPV has also been associated with more severe illness and lower respiratory tract symptoms such as pneumonia and bronchitis.


[Test principle]

This kit uses the double antibody-sandwich method to detect antigens. When an appropriate amount of specimen is added to the specimen well(s) of the test device, the specimen will move forward along the test device. If the specimen contains an antigen, the antigen combines with the antibody immobilized on the conjugate pad, and the immune complex forms a sandwich complex with another coated antibody which was coated on the test line, a visible colored line will show up, which indicates that the antigen is positive. The test device also contains a quality control line, regardless of whether there is a test line, the red quality control line should appear. If the quality control line does not appear, it indicates that the test result is invalid and you need to do the test again.

[Limitations of inspection methods]

1. This product is used for qualitative testing only and cannot indicate the level of antigen in the specimen.
2. Negative results may occur if the antigen titre in the specimen falls below the minimum detection limit of this kit.
3. This test kit is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor's diagnosis shall prevail.
4. Diagnosis and treatment can not only rely on this test result. Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.
5. The accuracy of the test depends on the quality of the swab sample, false negative results may be given following poor sampling.


Positive



Negative result:

If only quality control line C is colored, test line A or test line B or test line COVID-19 or test line RSV or test line ADV or test line HMPV are colorless, it means that no antigen of corresponding pathogen is detected, and the result is negative.


Negative



Invalid result:

If the quality control line C is not observed, it will be invalid regardless of whether there is test line A or test line B or test line COVID-19 or test line RSV or test line ADV or test line HMPV, and the test shall be conducted again.

Invalid



- Any failure to respect the test procedure may negatively impact the performance of the test and/or invalidate the test result.
- If the result of the test is negative, yet clinical symptoms persist, it is advised that you carry out additional tests using other clinical methods. A negative result at no time rules out the presence of antigens of the virus in the sample, as they may be present but at a level inferior to the minimum detection level of the test, or if the sample has been collected incorrectly.
- This test is not a substitute for a medical consultation, or for the result of a biological analysis carried out in a medical analysis laboratory.

[Warnings and Precautions]

- For *in vitro* diagnostic use.
- Read the instructions carefully before using the kit, and strictly control the reaction time. If you do not follow the instructions, you will get inaccurate results.
- Guard against moisture, do not open the aluminum platinum bag before it is ready for testing. Do not use it if the aluminum foil bag is damaged or the test device is damp.
- Please use it within the validity period.
- Balance all reagents and specimens to room temperature (15 ~ 30 °C) before use.
- Do not replace the components in this kit with components in other kits.
- Do not dilute the specimen when testing, otherwise you may get inaccurate results.
- The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
- The test methods and results must be interpreted in strict accordance with this specification.
- Negative results may occur if the antigen titer in the specimen falls below the minimum detection limit of this kit.
- Do not eat, drink or smoke in the area where the specimens or kits are handled. Wash hands thoroughly after finishing the tests.
- Do not perform the test in a room with strong air flow, ie. an electric fan or strong air-conditioning.
- Do not mix and interchange different specimens.
- Any serious incident occurring during the use of the equipment should be reported to the manufacturer and to the competent authorities of the Member State in which the user and/or the patient is based.

[Storage conditions & period of validity]

- The kit should be stored at 2-30 °C until the expiry date printed on the sealed pouch.
- After the foil pouch is unsealed, the test device should be used as soon as possible within one hour.
- The test device should be kept away from direct sunlight, moisture and heat.
- Do not freeze the test kit.

[Sample Transport and Storage]

After Swab specimens were collected, swab can be stored in extraction reagent provided with the kit. Freshly collected specimens should be processed as soon as possible, but no later than one hour after specimen collection. Specimen collected may be stored at 2-8°C for no more than 7 days.

[Quality Control]

Program control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient volume of the specimen.

[Performance index]

1. Physical characters

1.1 Appearance

The test should be clean and complete, no burr, no damage and non-pollution. The shell of the test cassette should be flat, the upper and lower covers should be evenly closed, and there should be no obvious gap. The inner test strip should be firmly attached without wobble. The extraction reagent (0.24% NaH₂PO₄ • 2H₂O + 3.04% Na₂HPO₄ • 12H₂O + 0.4% Tween 20 + 0.1% PC-300) should be free of foreign matter.

1.2 Size

The size of the inner strip should not be less than 2.5mm.

1.3 Liquid migration speed

It should not be less than 10mm/min.

2. Sensitivity

2.1.Limit of detection (LOD)

Category		LOD
Influenza A	Liao Ning/1183/2007 (H1N1)	2.2 × 10 ³ TCID ₅₀ /mL
	A/Victoria/3/75	1.16 × 10 ² TCID ₅₀ /mL
	A/HongKong/8/68	2.58 × 10 ⁴ TCID ₅₀ /mL
	A/Vietnam/1194/2004(H5N1)	2.98 × 10 ³ TCID ₅₀ /mL
	A/waterfowl/Chenhu/367-1/2021	2.05 × 10 ⁴ TCID ₅₀ /mL
	A/Wisconsin/67/2022	1.65×10 ³ TCID ₅₀ /mL
	A/Norway/31694/2022	1.49×10 ³ TCID ₅₀ /mL
	A/California/122/2022	2.89×10 ³ TCID ₅₀ /mL
Influenza B	A/Thailand/8/2022	1.96×10 ³ TCID ₅₀ /mL
	B/Phuket/3073/2013	1.78 × 10 ³ TCID ₅₀ /mL
	B/HongKong/574/2019	1.82 × 10 ³ TCID ₅₀ /mL
	Virus B/Florida/78/2015	7.8 × 10 ³ TCID ₅₀ /mL
	Jiang Xi/32/2000	2.9 × 10 ² TCID ₅₀ /mL
	B/1704	6.8 × 10 ² TCID ₅₀ /mL
	B/Brisbane/46/2015	2.05×10 ³ TCID ₅₀ /mL
	B/Austria/1359417/2021	2.67×10 ³ TCID ₅₀ /mL
B/Brisbane/1/2007	2.43×10 ³ TCID ₅₀ /mL	
RSV	RSV	10 ng/mL
RSV A	A2,ATCC VR-26	3.35 × 10 ³ TCID ₅₀ /mL
RSV B	WYZ9320,CH93(18)-18	3.75 × 10 ² TCID ₅₀ /mL
COVID-19		80 TCID ₅₀ /mL
ADV		10 ng/mL
HMPV		1 x10 ² TCID ₅₀ /mL

The LOD established with 1st International Standard for SARS-CoV-2 antigen NIBSC code: 21/368 (Version 1.0, Dated 01/11/2022) is 10 IU/ mL.

3. Specificity

3.1 Cross reaction

The cross-reactivity of this reagent is evaluated by testing a group of related pathogens, high-prevalence disease pathogens, and normal or pathogenic flora. The results prove that the product has no cross-reactivity, including Staphylococcus aureus, Streptococcus pneumoniae, Measles virus, Mumps virus, Bordetella pertussis, EB virus, Rhinovirus, Chlamydia pneumoniae, Human coronavirus OC43, Human coronavirus 229E, Human coronavirus NL63, MERS coronavirus, Legionella pneumophila, Parainfluenza virus, Mycobacterium tuberculosis, Mycoplasma pneumoniae, Enterovirus, Staphylococcus epidermidis, Streptococcus pyogenes, Haemophilus influenzae, Candida albicans, Pseudomonas aeruginosa, Streptococcus salivarius. Moreover, there is no cross-reactivity among the six viruses tested by this kit.

3.2 Interfering substances

Common interfering substances in the sample (such as blood, mucus, and pus) as well as common medications for respiratory diseases and nasal sprays do not affect the test results.

4. Clinical performance

4.1 COVID-19 Rapid Test:

A total of 1011 samples were collected in this study, of which 103 were positive and 908 were negative. The statistics of the study results are shown in the following table:

Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE COVID-19 Test		POS	NEG	Total	PPA	98.06%	93.19%	99.47%
	POS	101	0	101	NPA	>99.9%	99.58%	100%

	NEG	2	908	910	Total compliance rate	99.80%
	TOTAL	103	908	1011		

Diagnostic sensitivity: 98.06%; Diagnostic specificity: >99.9%; Total compliance rate: 99.80%

4.2 Influenza A Rapid Test:

A total of 1011 samples were collected in this study, of which 107 were positive and 904 were negative. The statistics of the study results are shown in the following table:

Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE Influenza A Test		POS	NEG	Total	PPA	98.13%	93.44%	99.49%
	POS	105	0	105	NPA	>99.9%	99.58%	100%
	NEG	2	904	906	Total compliance rate		99.80%	
	TOTAL	107	904	1011				

Diagnostic sensitivity: 98.13%; Diagnostic specificity: >99.9%; Total compliance rate: 99.80%

4.3 Influenza B Rapid Test:

A total of 1011 samples were collected in this study, of which 101 were positive and 910 were negative. The statistics of the study results are shown in the following table:

Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE Influenza B Test		POS	NEG	Total	PPA	95.05%	88.93%	97.87%
	POS	96	0	96	NPA	>99.9%	99.58%	100%
	NEG	5	910	915	Total compliance rate		99.51%	
	TOTAL	101	910	1011				

Diagnostic sensitivity: 95.05%; Diagnostic specificity: >99.9%; Total compliance rate: 99.51%

4.4 RSV Rapid Test:

A total of 1011 samples were collected in this study, of which 103 were positive and 908 were negative. The statistics of the study results are shown in the following table:

Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE RSV Test		POS	NEG	Total	PPA	95.15%	89.14%	97.91%
	POS	98	0	98	NPA	>99.9%	99.58%	100%
	NEG	5	908	913	Total compliance rate		99.51%	
	TOTAL	103	908	1011				

Diagnostic sensitivity: 95.15%; Diagnostic specificity: >99.9%; Total compliance rate: 99.51%

4.5 ADV Rapid Test:

A total of 1011 samples were collected in this study, of which 101 were positive and 910 were negative. The statistics of the study results are shown in the following table:

Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE ADV Test		POS	NEG	Total	PPA	95.05%	88.93%	97.87%
	POS	96	0	96	NPA	>99.9%	99.58%	100%
	NEG	5	910	915	Total compliance rate		99.51%	

	TOTAL	101	910	1011		
--	-------	-----	-----	------	--	--

Diagnostic sensitivity: 95.05%; Diagnostic specificity: >99.9%; Total compliance rate: 99.51%
















4.6 Human metapneumovirus (HMPV) Rapid Test:

A total of 1011 samples were collected in this study, of which 100 were positive and 911 were negative. The statistics of the study results are shown in the following table:


Reference RT-PCR analysis							95% Wilson Score CI	
							LCI	UCI
DEEPBLUE HMPV Test		POS	NEG	Total	PPA	97.00%	91.55%	98.97%
	POS	97	0	97	NPA	>99.9%	99.58%	100%
	NEG	3	911	914	Total compliance rate		99.70%	
	TOTAL	100	911	1011				

Diagnostic sensitivity: 97.00%; Diagnostic specificity: >99.9%; Total compliance rate: 99.70%

[Index of Symbols]

	In vitro diagnostic medical device		Do not re-use		Keep away from sunlight
	Use-by date		Consult instructions for use or consult electronic instructions for use		Date of manufacture
	Caution		Manufacturer		Do not use if package is damaged and consult instructions for use
	Temperature limit		Batch code		Contains sufficient for <n> tests
	Authorized representative in the European Union		Keep dry		CE Mark

 ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO., LTD.
No. 777 Jimingshan Road, High-Tech Development Zone, 230088 Hefei, Anhui, PEOPLE'S REPUBLIC OF CHINA
E-mail:sales@dblumedical.com

 Mega Eurostar Sp. z o. o.
Obrzeźna, 5 lok. XIP/1 02-691 Warsaw Poland

The Safety Data Sheet (SDS) can be obtained via e-mail.

Specification	REF
1 piece per box	CFRAH1NST-1
2 pieces per box	CFRAH1NST-2
5 pieces per box	CFRAH1NST-5
10 pieces per box	CFRAH1NST-10



Scan QR code for IFU in different languages.

